

Unter Haßlen 20 · D-78532 Tuttlingen
Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50
e-mail: sales@rz-medizintechnik.com
<http://www.rz-medizintechnik.com>



UROLOGY Disposables



10/09

| RZ Hydrophilic Guide Wire Nitinol

Advantages:

- Nitinol Core
- Extruded PU Coating
- Hydrophilic Surface
- Excellent Visibility on X-Rays



| 253-832-150

0.032"

straight tip, length: 150 cm
one side rigid, one side flexible

| 253-832-151

0.032"

curved tip, length: 150 cm
one side rigid, one side flexible

| 253-835-150


0.035"

straight tip, length: 150 cm
one side rigid, one side flexible

| 253-835-151

0.035"


curved tip, length: 150 cm
one side rigid, one side flexible



| RZ Stone Baskets

Advantages:

- optimal shaping results in high opening force
- various designs
- different tips and lengths *
- suitable for 1.9 – 7.0 Charrier tubes



| 253-960-025


Stone Basket, 2.5 Charr., sterile for
single use, package of 5 pieces

* RZ Medizintechnik offers a wide range of different Stone Basket Instruments. Feel free to ask for your desired version.

RZ Ureteral Balloon Dilator

Advantages:

- non compliant
- short flexible, radiopaque tip
- puncture resistant balloon material
- high rated burst pressure
- low profile
- fast deflation
- new flexible catheter shaft
- clear catheter outer shaft
- flexible balloon



RZ Reference No.	Balloon Ø	Balloon Length	Catheter Ø	Rated Burst Pressure	Guide Wire
Shaft Length: 75 mm	mm / CH	mm	mm / CH	atm (RBP)	Inch
254-904-040	4.0 / 12	40	2.0 / 6.0	20	0.035
254-904-060	4.0 / 12	60			
254-904-100	4.0 / 12	100			
254-905-040	5.0 / 16	40			
254-905-060	5.0 / 16	60			
254-905-100	5.0 / 16	100			
254-906-040	6.0 / 18	40			
254-906-060	6.0 / 18	60			
254-906-100	6.0 / 18	100			
254-907-040	7.0 / 21	40			
254-907-060	7.0 / 21	60			
254-907-100	7.0 / 21	100			
254-908-040	8.0 / 24	40			
254-908-080	8.0 / 24	80			
254-910-040	10.0 / 30	40			
254-910-060	10.0 / 30	60			
254-910-080	10.0 / 30	80			

RZ percutaneous dilators for renal aditus

RZ reference: 253-900-170

sterile



DILATORS



SHEATHS

Description	Length (mm)	Fr
Facial Dilator 6 Fr	300	6
Facial Dilator 8 Fr	300	8
Facial Dilator 10 Fr	300	10
Facial Dilator 12 Fr	300	12
Facial Dilator 14 Fr	300	14
Facial Dilator 16 Fr	300	16
Facial Dilator 18 Fr	300	18
Renal Dilator 20 Fr	300	20
Renal Dilator 22 Fr	300	22
Renal Dilator 24 Fr	300	24
Renal Dilator 26 Fr	300	26
Renal Dilator 28 Fr	300	28
Renal Dilator 30 Fr	300	30
Introducer Catheter	700	7

Fr (ID)	Fr (OD)	Length (mm)
26	30	180
28	32	180
30	34	180

Available soon: Sterile RZ Nephrostomy Set, Pigtail & Malecot