

Unter Haßlen 20 · D-78532 Tuttlingen  
Phone: +49 (0)7462/9470-0  
Fax: +49 (0)7462/9470-50  
e-mail: sales@rz-medizintechnik.com  
<http://www.rz-medizintechnik.com>



## SPINE PREP Preparation Set



10.12

# SPINE PREP DISC PREPARATION SET

The new RZ SPINE PREP instruments facilitate discectomy for anterior lumbar procedures. This comprehensive set allows access to posterolateral corners by its optimized working length which is also suitable for large body masses. Maximum control is assured by the handles with double handed grip.



| 259-156-010  
Box Curette, 6mm x 10mm, teflon handle, 430mm length



| 259-156-020  
Endplate Elevator, 20mm width, teflon handle, 430mm length



| 259-156-001 3.5mm x 4.5mm  
| 259-156-002 5.5mm x 8.5mm  
| 259-156-003 7.5mm x 11.5mm  
Bone Curette, angled, oval, teflon handle, 430mm length



| 259-156-025  
Ball Tip Probe 90°, 300mm length





- | 250-502-330 2mm
- | 250-504-330 4mm
- | 250-506-330 6mm

**Ferris-Smith- Kerrison**  
Laminectomy Punch, 40° upwards cutting, 330mm



- | 250-303-330 3 x 10mm
- | 250-304-330 4 x 10mm
- | 250-306-330 6 x 10mm

**Ferris-Smith**  
Rongeur for Thoracic Spine Surgery, upwards cutting, 330mm working length





| 250-801-360  
SYPERT Bone Rongeur,  
straight, smooth jaws  
360mm



| 250-802-360  
SYPERT Bone Rongeur,  
slightly curved, smooth jaws  
360mm



| 250-803-360  
SYPERT Bone Rongeur,  
straight, serrated jaws  
360mm