



progress in surgery

GYNECOLOGY
LEEP / LLETZ
PROGRAM



GYNECOLOGY

GYNECOLOGY LEEP / LLETZ PROGRAM

With this range of blue coated instruments and complimentary devices, RZ provides all necessary products for the LEEP (Loop Electrical Excision Procedure) and LLETZ (Large Loop Excision of the Transformation Zone) operation methods. The surgeon can choose between various models of the vaginal specula, forceps, hooks and single use as well as reusable electrodes. A complete Set of G400 with its accessories and a smoke evacuation unit round the set off - ready to operate.

These procedures are done to remove abnormal, precancerous or cancerous cells from the cervix or the endocervix, to diagnose cancer, to try to answer questions not answered by a colposcopy and to better understand the results of a PAP test in general.

Benefits of these methods are the finding and removing of precancerous cells to prevent cancer, the determination and diagnosis of the patients problem in general as well as the possibility to perform the procedure in the doctors office because it does not require a general anesthetic, only a local or regional one.

The speculas are equipped with a built-in smoke evacuation port to remove the upcoming smoke completely and provide a smoke-free operating field. With the durable insulating blue coating of all instruments we ensure the highest level of protection for the patient as well as for the operating team from burns or shockings. The coated instruments can be sterilized using the common cleaning cycles.

The four-sided specula come up with two additional advantages: an even better view and access to the cervix for the surgeon and a higher safety level for the patient concerning vaginal wall burns.

The four-sided specula combines the two-sided specula and the lateral retractor.



GYNECOLOGY



Graves

350-131-001 Small

350-131-002 Medium

350-131-003 Large



Pederson

350-141-002 Medium

350-141-003 Large



Collin

350-111-003 Medium

350-111-004 Large



Weisman-Pederson

350-140-000 Right, medium

350-141-000 Left, medium



Weisman-Graves

350-138-000 Left, medium

350-142-002 Left, large

350-137-000 Right, medium

350-142-001 Right, large

GYNECOLOGY LEEP / LLETZ PROGRAM



350-510-070




350-514-070



350-518-070

VuMar

350-510-070	Medium
350-512-070	Large
350-514-070	XL, Graves medium
350-516-070	XL, Graves large
350-518-070	XL, Pederson medium
350-520-070	XL, Pederson large



Clinton

350-133-001	Wide (Graves)
350-133-002	Narrow (Pederson)



350-790-001



350-790-003



350-790-004

Side Wall Retractor

350-790-001	4-sided Pederson speculum, medium
350-790-002	4-sided Pederson speculum, large
350-790-003	Redundant wall retractor
350-790-004	Lateral Wall Retractor, closed shanks
350-790-005	Lateral Wall Retractor, open shanks



Kogan

350-780-280 Ratchet, standard tips

350-782-280 Ratchet, narrow tips

350-784-280 Set screw, standard tips


350-786-280 Set screw, narrow tips



Schroeder

150-150-250 Ring forceps

350-401-250 Single tooth, coated tips



Hooks

350-768-240 Two-prong hook

350-769-240 Three-prong hook

350-767-240 Emmet hook, fully curved

350-767-241 Emmet hook, angled

350-767-242 Emmet hook, double angled

350-770-250 Iris hook

GYNECOLOGY LEEP / LLETZ PROGRAM



Serrated, elevated shanks

Potts-Smith dressing forceps

130-738-181	180mm, 7"
-------------	-----------

130-738-211	210mm, 8 ¹ / ₄ "
-------------	--

130-738-251	250mm, 10"
-------------	------------

130-738-301	300mm, 12"
-------------	------------



1x2 teeth, elevated shanks

Potts-Smith tissue forceps

130-741-181	180mm, 7"
-------------	-----------

130-741-211	210mm, 8 ¹ / ₄ "
-------------	--

130-741-251	250mm, 10"
-------------	------------

130-741-301	300mm, 12"
-------------	------------



Design your individual forceps on

www.dreamtweezer.com

The innovative Forceps Configurator for special requirements or special anatomical conditions.



Curved Tip

AxB	single use	reusable
5 x 5	-	700-750-055
10 x 5	-	700-750-105
10 x 10	700-755-105	700-750-100
15 x 8	-	700-750-158
15 x 10	700-755-158	-
15 x 12	700-755-152	700-750-152
15 x 15	700-755-155	700-750-155
20 x 8	-	700-750-208
20 x 10	700-755-210	-
20 x 12	700-755-212	-
20 x 15	700-755-215	700-750-215
20 x 18	700-755-251	-
20 x 20	700-755-208	700-750-220
25 x 10	-	700-750-251
25 x 22	-	700-750-225



Square Tip

AxB	single use	reusable
10 x 4	-	700-750-104
10 x 8	-	700-750-108



Electrode Sterilizing Case

400-607-203 203 x 114 x 30 mm

Electrodes

	single use	reusable
Needle Electrode \varnothing 0.5mm	-	700-760-510
Needle Electrode \varnothing 0.7mm	-	700-760-520
Ball Electrode 3mm / 2mm	700-755-003	700-760-002
Ball Electrode 5mm / 4mm	700-755-005	700-760-004
Electrode Extension (140mm)	-	700-760-900
Electrode Extension (90mm)	-	700-760-910



Needle Electrode



Ball Electrode



Electrode Extension

GYNECOLOGY LEEP / LLETZ PROGRAM



700-350-400 230V
700-350-410 110V
Electrosurgical Generator G400

- Detects and identifies connected instrument automatically
- Powerstart System simplifies the beginning of cutting process
- Automatic power regulation
- 105 different memory locations in 7 groups
- Various output configuration of the sockets available
- Footswitch (optional)
- Available with Argon and / or the RZeal



700-350-405
2-Button-Footswitch, length: 5m

700-350-402
*2-Button-Footswitch and 1 Vario-Switch-Button, length: 5m



700-503-003
Cable for Disposable Neutral Electrode
5m, only RZ compatible



700-503-000
Disposable Neutral Electrode for Adults,
hydrogel, split, for adults or children
10x5 pcs. per package



700-700-020
Hand Control Handle for Electrodes,
with 2 buttons, cable length: 5m




700-700-018
Footswitch Handle for Electrodes,
cable length: 5m





166-957-000
ATMOS Smoke evacuation basic unit
incl. Footswitch and connecting tube
for the olive on the specula (2.3m)



GYNECOLOGY LEEP / LLETZ PROGRAM



GYNECOLOGY





progress in surgery

RZ Medizintechnik GmbH

Unter Hasslen 20
D-78532 Tuttlingen

Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50

sales@rz-medizintechnik.com
www.rz-medizintechnik.com

Media Code Number: 021580