

Unter Haßlen 20 · D-78532 Tuttlingen
Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50
e-mail: sales@rz-medizintechnik.com
<http://www.rz-medizintechnik.com>




BiPOstyle Bipolar Scissors and Forceps



13/09

BIPOLAR MICRO SCISSORS





RZ Micro Scissors are especially made for precise cutting when it comes to highly delicate surgery. Its thin scissors blades are a major advantage e.g. for plastic and reconstructive surgery.



| 120-742-130
Bipolar Scissors 13.0cm - 5", curved blades, slim and thin scissors blades, for highly precise working




| 120-742-140
Bipolar Scissors 14.0cm - 5¹/₂", curved blades, slim and thin scissors blades, for highly precise working




| 120-742-175
Bipolar Scissors 17.5cm, - 7", curved blades, slim and thin scissors blades, for highly precise working


BIPOLAR SCISSORS AND CABLES



- | 120-740-180 RZ Bipolar Scissors, curved blades, 18.0cm, 7"
- | 120-740-210 RZ Bipolar Scissors, curved blades, 21.0cm, 8 1/4"
- | 120-740-230 RZ Bipolar Scissors, curved blades, 23.0cm, 9"
- | 120-740-280 RZ Bipolar Scissors, curved blades, 28.0cm, 11"



- | 120-741-180 RZ Bipolar Scissors, curved slim blades, 18.0cm, 7"
- | 120-741-210 RZ Bipolar Scissors, curved slim blades, 21.0cm, 8 1/4"
- | 120-741-230 RZ Bipolar Scissors, curved slim blades, 23.0cm, 9"
- | 120-741-280 RZ Bipolar Scissors, curved slim blades, 28.0cm, 11"



Bipolar Cables for use with RZ Bipolar Scissors



Valleylab



Martin



Erbe



Banana plugs



3m	700-736-003	700-737-003	700-738-003	700-739-003
5m	700-736-005	700-737-005	700-738-005	700-739-005

RZ BIPOLAR CLAMPS



I 700-730-230 RZ Bipolar Clamp, curved jaws,
23.0 cm, 9"



I 700-735-230 RZ Bipolar Clamp, curved jaws, with ratchet,
21.0 cm, 8¹/₄"



700-731-500 Bipolar Cable for RZ Bipolar Clamp,
with 2 banana plugs, 5m length



The RZ bipolar scissors-forceps enables powerful grasping and compression of tissue, combined with bipolar coagulation or sealing of blood vessels. Now, tissue can be cut without changing the instrument.



I 700-120-190 RZ bipolar scissors-forceps, straight, with serrated
jaws, without ratchet, 19.0 cm, 7¹/₂"

I 700-120-140 RZ bipolar scissors-forceps, straight, with serrated
jaws, without ratchet, 14.0 cm, 5¹/₂"



I 700-121-190 RZ bipolar scissors-forceps, **curved**, with serrated
jaws, without ratchet, 19.0 cm, 7¹/₂"

I 700-121-140 RZ bipolar scissors-forceps, **curved**, with serrated
jaws, without ratchet, 14.0 cm, 5¹/₂"