

Unter Haßlen 20 · D-78532 Tuttlingen  
Phone: +49 (0)7462/9470-0  
Fax: +49 (0)7462/9470-50  
e-mail: sales@rz-medizintechnik.com  
http://www.rz-medizintechnik.com



## CASPAR MICRO Lumbar Discectomy Retractor



03/11

## CASPAR MICRO LUMBAR DISCECTOMY RETRACTOR

- Specially designed Retractor Set allows a micro surgery approach to Lumbar Disc
- Jointed arms folding along the patients back allow unhindered access around the wound
- All speculas and blades have a black surface finish to prevent glare
- Fenestrated blades on the specula reduce the risk of movement during surgery



| 259-185-100  
Caspas Counter Retractor, with black finished hinged arms



| 259-185-101  
Blade Ejector Forceps



| Caspar Specula with black finished surface

	blade length
259-185-040	40mm
259-185-045	45mm
259-185-050	50mm
259-185-055	55mm
259-185-060	60mm
259-185-065	65mm
259-185-070	70mm



Fast, easy and safe locking of different types of blades with SnapLock mechanism

# LATERAL AND MEDIAL BLADES

## Lateral Blade with SnapLock, black finish

	blade length
259-185-235	35mm
259-185-240	40mm
259-185-245	45mm
259-185-250	50mm
259-185-255	55mm
259-185-260	60mm
259-185-265	65mm
259-185-270	70mm
259-185-275	75mm
259-185-280	80mm
259-185-285	85mm



## Medial Blade, black finish

	blade length
259-185-340	40mm
259-185-345	45mm
259-185-350	50mm
259-185-355	55mm
259-185-360	60mm
259-185-365	65mm
259-185-370	70mm



Please ask for our exclusive line of spinal Instruments